


# **AcryRock**

*We take care of your smile*



DENTI ARTIFICIALI IN RESINA - ACRYLIC TEETH



## Scala Colori - Shade Guide

Acry Rock V
A1 - A2 - A3 - A3,5 - A4
B1 - B2 - B3 - B4
C1 - C2 - C3 - C4
D2 - D3 - D4
E1
F1 - F2





**Tabella di Articolazione**  
**Articulation table**

S 13	I 13	D 35
S 14	I 13	D 35
S 15	I 38	D 37 - D 43
S 17	I 41	D 39 - D 44
S 18	I 18	D 39 - D 44
S 25	I 18	D 39
S 26	I 46	D 39 - D 45
S 32	I 41 - I 42	D 39 - D 44
S 33	I 46	D 45
S 37	I 37	D 36 - D 42
S 38	I 38	D 37 - D 43
S 39	I 51	D 35
S 40	I 40	D 36 - D 42
S 41	I 41	D 39 - D 44
S 42	I 42	D 39 - D 44
S 47	I 18 - I 53	D 39
S 48	I 46	D 45
S 51	I 51	D 35
S 52	I 52	D 36 - D 42
S 53	I 53	D 39 - D 44
S 54	I 53	D 39 - D 44
S 55	I 66	D 45 - D 46
S 56	I 46	D 45
S 60	I 60	D 33
S 61	I 61	D 36 - D 42
S 62	I 61	D 36 - D 42
S 63	I 63	D 39 - D 44
S 64	I 64	D 41
S 65	I 65	D 41
S 66	I 66	D 44
S 67	I 67	D 45
S 68	I 68	D 46





# AcryRock

*We take care of your smile*




Codice Prodotto Product Code	
---------------------------------	--

0010001	ACRY ROCK 6 UPPER
0010002	ACRY ROCK 6 LOWER
0010003	ACRY ROCK 8 UPPER
0010004	ACRY ROCK 8 LOWER
0020001	ACRY ROCK SET 28
0020002	LIVE MOULD CHART ACRY ROCK
0030001	SHADE GUIDE ACRY ROCK

0010013	ACRY ROCK V. 6 UPPER
0010014	ACRY ROCK V. 6 LOWER
0010015	ACRY ROCK V. 8 UPPER
0010016	ACRY ROCK V. 8 LOWER
0020003	ACRY ROCK V. SET 28
0020004	LIVE MOULD CHART ACRY ROCK V.
0030002	SHADE GUIDE ACRY ROCK V.



Acry Rock comply with ISO 22112:2005 - Type 1 - Type 2



**A:** Anatomico - *Anatomic*  
**G:** Geometrico - *Geometric*

**D33** (1) 26,9 (2) 30,5 (3) 7,2 (4) 8,1



**D35** (1) 30,7 (2) 33,0 (3) 8,5 (4) 5,5



**D42** (1) 28,7 (2) 30,6 (3) 8,2 (4) 7,8



**D36** (1) 28,8 (2) 31,1 (3) 8,2 (4) 8,6



**D43** (1) 29,6 (2) 31,4 (3) 8,3 (4) 8,1



**D37** (1) 30,7 (2) 33,2 (3) 8,4 (4) 8,8



**D44** (1) 31,1 (2) 33,0 (3) 8,5 (4) 8,6



**D39** (1) 30,1 (2) 35,0 (3) 8,5 (4) 8,5



**D45** (1) 31,0 (2) 33,3 (3) 8,8 (4) 8,6



**D41** (1) 30,0 (2) 33,5 (3) 8,6 (4) 8,9



**D46** (1) 33,6 (2) 35,6 (3) 9,6 (4) 9,2





Dental Manufacturing S.p.A.  
Via Ca' Mignola Nuova, 1699  
45021 Badia Polesine (RO) - Italy  
Tel. +39.0425.51628 - +39.0425.590587  
Fax +39.0425.590156  
[www.ruthinium.it](http://www.ruthinium.it) [info@ruthinium.it](mailto:info@ruthinium.it)